

---

## **ACEI Antibody**

**Subcategory:** Rabbit Polyclonal Antibody

**Cat. No.:** 250451

**Unit:** 0.1 mg

**Description:**

Angiotensin I converting Enzyme Inhibitor (ACEI)

**Isotype:** Rabbit Ig

**Applications:** E, WB, IHC

**Species Reactivity:** H, M, R, D

**Format:** Each vial contains 0.1 mg IgG in 0.1 ml (1 mg/ml) of PBS pH7.4 with 0.09% sodium azide. Antibody was purified by Protein-G affinity chromatography.

**Alternate Names:** Tuna muscle-derived angiotensin-converting enzyme inhibitory peptide; ACEI; PTHIKWGD; angiotensin-converting enzyme inhibitor peptide; Angiotensin-converting enzyme inhibitor

**Accession No.:** P18691

**Antigen:** KLH-conjugated tuna ACEI.

**Application Notes:** E: 1:500-1:1,000; WB: 1:200-1:500; IHC: 1:200-1:500

**Storage:** Store at -20°C. Minimize freeze-thaw cycles.

Product is guaranteed one year from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.